



The right probe for any application

The ACCUVIX V10 optional probes provide exceptional imaging performance in all applications. Advanced technology of automatic 3D volume reconstructing and expanded functional bandwidth delivers increased amount of clinical information and qualified diagnosis.



www.medison.com | info@medison.com



see it all
MEDISON

PG-MD-V10-ICM-09.01.2008



























Probe Guide

Version 1.0

see it all
MEDISON

Probe Guide

																								
TRANSDUCER	C2-5EL	C2-6IC	C3-7IM	EC4-9IS	NEV4-9ES	NER4-9ES	EV4-9ED	ER4-9ED	L6-12IS	L5-13IS	L8-15IS	L7-16IS	L5-12/50EP	L4-7EL	C4-9/10ED-N	3D2-6ET	3D4-8ET	EC5-9EK	P2-4AC	P3-5AC	P2-4BA	P3-8CA	2.0CW	4.0CW
Array Type	Curved	Curved	Curved	Curved	Curved	Curved	Curved	Curved	Linear	Linear	Linear	Linear	Linear	Linear	Curved	Curved	Curved	Curved	Phased	Phased	Phased	Phased	CW	CW
Sector Angle / Field of View	76°	58°	78°	149°	150°	150°	148°	148°	38mm	38mm	38mm	38mm	50mm	38mm	150°	84°	84°	128°	19mm	14mm	80mm	15mm	N/A	N/A
Radius of Curvature	40mm	60mm	50mm	10mm	10mm	10mm	10mm	10mm	Flat	Flat	Flat	Flat	Flat	Flat	10mm	40mm	40mm	10mm	Flat	Flat	Flat	Flat	N/A	N/A
Broadband Frequency Range	2MHz~5MHz	2MHz~6MHz	3MHz~7MHz	4MHz~9MHz	4MHz~9MHz	4MHz~9MHz	4MHz~9MHz	4MHz~9MHz	6MHz~12MHz	6MHz~12MHz	8MHz~15MHz	8MHz~15MHz	5MHz~12MHz	4MHz~7MHz	4MHz~9MHz	2MHz~6MHz	4MHz~8MHz	5MHz~9MHz	2MHz~4MHz	3MHz~5MHz	2MHz~4MHz	3MHz~8MHz	2MHz	4MHz
APPLICATION																								
Abdomen	●	●	●													●	●							
Gynecology	●	●	●	●	●	●	●	●								●	●	●						
OB	●	●	●	●	●	●	●	●								●	●	●						
OB Early	●	●	●	●	●	●	●	●								●	●							
Fetal Heart	●	●	●	●	●	●	●	●								●	●							
Musculoskeletal									●	●	●	●	●	●										
Urology	●	●	●	●	●	●	●	●										●						
Pediatric									●	●	●	●	●	●	●				●	●	●	●		
Small Part									●	●	●	●	●	●										
Breast									●	●	●	●	●	●										
Vascular									●	●	●	●	●	●	●									
Cardiology																			●	●	●	●	●	●
TCD																			●	●	●	●	●	●
Neonatal															●				●	●	●	●	●	●
Contrast Agent																								
NEEDLE GUIDE	DISPOSABLE	DISPOSABLE	REUSABLE	REUSABLE	REUSABLE/ DISPOSABLE	REUSABLE/ DISPOSABLE	N/A	N/A	DISPOSABLE	DISPOSABLE	DISPOSABLE	DISPOSABLE	DISPOSABLE	DISPOSABLE	N/A	DISPOSABLE	DISPOSABLE	REUSABLE	N/A	N/A	N/A	N/A	N/A	N/A
Biopsy Kit	BP-KIT-010	BP-KIT-009	BP-C3-7IM	BP-KIT-003	REUSABLE BP-KIT-019 DISPOSABLE BP-KIT-018	REUSABLE BP-KIT-019 DISPOSABLE BP-KIT-018			BP-KIT-012	BP-KIT-012	BP-KIT-012	BP-KIT-012	BP-KIT-020	BP-KIT-014		BP-KIT-017	BP-KIT-017	VF-BP-PEC47						
	